Pt

**Unimatrix**

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### Unimatrix

The Unimatrix system is indicated for Class II restorations, and composed by preformed sectional matrices and Unimatrix rings. The main advantages of this system offers are the easiness in the reconstruction of contact points and natural tooth anatomy. Unimatrix rings are equipped with silicone tips for a better adaptation and to bring more comfort to the patient. **Use:** Select, among the 3 matrix sizes, the one that best adapts the cavity to be restored. Install the matrix in the proximal surface of the tooth to be restored. When it concerns to a MOD cavity, it is advisable to restore one surface at a time. Install the wedge. TDV elastic wedges improve the adaptation of the matrix; in this case, silicone tips are not used. Install the Unimatrix ring. The use or not of the ring and silicone tips, as well as setting up manner depends of the personal preferences and the clinical situation. Concluding the restoration, remove first the Unimatrix ring, then the wedges and lastly the matrix. To obtain consistent contact points, it is recommended to separate the adjacent teeth by using wedges or Unimatrix ring, from the anesthesia point until the installation of the matrix, as well as, the use of the Contact + instrument during the composite polymerization. **Sizes:** Available in small, medium and large sizes.

Ref.	Code	Content / Contenu
4122	0003020006	Mini-Kit - 25 assorted matrices + 1 ring + 4 tubes / Mini-Kit - 25 matrices assorties + 1 ring + 4 protecteurs
4124	0003020005	Kit - 50 assorted matrices + 2 rings + 8 tubes / Kit - 50 matrices assorties + 2 rings + 8 protecteurs
4125	0003020015	Refill - 50 assorted matrices / Recharge - 50 matrices assorties
4125G	0003020016	Refill - 50 large matrices / Recharge - 50 matrices grandes
4125M	0003020017	Refill - 50 medium matrices / Recharge - 50 matrices moyennes
4125P	0003020018	Refill - 50 small matrices / Recharge - 50 matrices petites
4126	0003020019	Unimatrix Rings - 2 pieces + 8 tubes / Rings Unimatrix - 2 unités + 8 protecteurs
4120	0003020022	Kit - 100 assorted matrices + 2 rings + 8 tubes / Kit - 100 matrices assorties + 2 rings + 8 protecteurs
4121	0003020029	Refill - 100 assorted matrices / Recharge - 100 matrices assorties
4121G	0003020030	Refill - 100 large matrices / Recharge - 100 matrices grandes
4121M	0003020031	Refill - 100 medium matrices / Recharge - 100 matrices moyennes
4121P	0003020032	Refill - 100 small matrices / Recharge - 100 matrices petites
4159	0003020063	Refill - Silicone tubes 20 units / Recharge - Protecteurs de silicone 20 unités
4160	0003020064	Refill - Silicone tubes 50 units / Recharge - Protecteurs de silicone 50 unités

**Warnings, precautions and restrictions:** Professional use. The matrices are for single use only. Its reuse can cause cross contamination. When there is proven allergy to any one of its components, the product must not be used. It is advisable to protect the papilla by any method to avoid accidents. It is advisable that the use of the product is always accompanied by absolute isolation of the operative field, as objects can be swallowed or breathed in by accident. The rings and the silicone tips must be sterilized separately. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. **Cleaning and sterilization:** The rings and the silicone tips are sterilizable in autoclave (126°C - 16min / 134°C - 18min). Before sterilization, proceed with the craning using water and enzymatic detergent. Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carbohydrase) and with pure pH range between 6.5 - 7.5. Prepare the solution according to the manufacturer's instructions. The water should be potable. For a better effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the products and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the pieces with absorbent paper or disposable wipes. Inspection: Check the product before each use, observing whether there are cracks or breakage, and do not use in these cases. These conditions are indicative of fragility of the material, whose strength can be reduced during clinical use. **Composition:** Matrices: Stainless steel. Ring: Nickel plated carbon steel. Tubes: Silicone. **Storage:** To be stored in a cool and dry place, protected from sunlight.

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Obelis S.A.