

Prophylactic Strip is an indispensable aid in can cause injuries to patient, we recommend prophylaxis procedures, since it easily accesses the interproximal region, where its active face effectively affect the product's performance. Any serious provides the removal of dental calculus. The circular incident related to the device should be micro-blades, located around the holes in one face reported to the manufacturer and the of the strip, has as main characteristic the action against calculus, thus preserving the enamel and restorative materials integrity

Indications: Prophylactic Strips is indicated for clinical use, removing dental plague and tartar. Use: Insert the strip in the interproximal space,

placing its active face turned to the desired tooth, Hold the strip by the ends, pull it against the tooth and apply back and forth longitudinal movements until achieve the desired effect. Remove the strip. Dimensions: 100mm x 4mm.

Warnings, precautions and restrictions: Professional use only. Single use only, Its reuse may cause cross contamination. When there is proven allergy to any of its component, the product Resp. Técn: Gisele K, Gruber must not be used. The strip in contact with the gum

competent authority of the Member State in which the user and/or patient is established. Storage: To be stored in a dry place, fresh and sheltered from sunlight. Composition: Stainless steel.

Dispose in accordance with applicable local

Ref. | Code Content 0.3076 0009010101 10 strips of 100mm x 4mm







ECREP Obelis S.A. Brl. Général Wahis 53 1030 Brussels, Belgium Tel: +(32) 2.732.59.54 E-mail: mail@obelis.net

Manufacturer: TDV Dental Ltda. CNPJ 81,591,786/0001-60 XV de Novembro, 9944 | CEP 89107-000 Pomerode/SC - Brasil Reg. ANVISA 10291220068 t. 55 47 3395-6115 info@tdv.com.br | www.tdv.com.br

Versão 07 - 06/2021



interproximal tartar removal 10 Non-abrasive steel strip with circular micro blades for

Prophylactic Strips

