Resinlay Pattern Photo - Pasta / Paste

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Resinlay Pattern Photo - Pasta / Paste

Light cured universal modelling acrylic resin in paste.

Resinlay Pattern Photo is a light-curing acrylic resin developed for universal modelling and casting patterns. Available in gel and paste consistencies, both versions can be used in associated techniques, allowing a wide number of applications by a complete polymerization and burning clearly leaving no residues, even on thicker layers. Due to its consistency, the paste version facilitates the sculpting and modelling on techniques that requires application of larger quantities.

Indication: This product is indicated for use by both, dentist and prosthetic technician. Among its main indications are: **Dentist:** Fixation and union of transfers in cases of multiple prostheses on two or more implants; Copings pattern sculpture; Caps for modeling; Fixation of temporary splint for soldering on bars for protocols or prosthetic structures of fixed prosthesis; Bite registration; Lucia jig; Post and Core impression by direct technique (associated use to gel); Fabrication and adjustment of guides: surgical, for positioning and insertion of provisional crowns; Reliefs. **Prosthesis Technician:** Preparation of body test for: implant bar, veneers, prostheses structure, copings; Ceramic structure; Abutments; Post and Core

Impression (associated use to gel); Copings for modeling; Reliefs.

Use: Before application, isolate the area with an insulating liquid (when used in plaster models) or Vaseline (when applied intraorally). This procedure allows the creation of a thin layer on the model, preventing the resin from sticking to the model while the impression is taken. It also makes easier the removal of the light-cured dental piece. And, it is still recommended to use solid Vaseline on hands and instruments, considering that the material has some adhesion. Open the cap and take the amount required for the procedure. Then, close the syringe, once the material exposure to the light for extended periods can compromise the product's functionality and performance. Before application, manipulate the product with your hands for a few seconds, to make it the proper consistency for application. Proceed with the product application. It is suggested to the applications be performed by layers up to 2 mm thick. In case of construction of a thicker piece, it is recommended to apply overlapping intermediate layers. Before forming the next layer, light cure for approximately 20 seconds without removing the dispersion layer. Emphasizing that the material's curing can be verified by changing the color from green to blue. For implant transfer copings, it is recommended to cut the bar between the posts and repair with Resinlay Pattern Photo - Gel, to minimize the risks of material contraction. Perform the finishing of the dental model, when applicable. The dispersion layer can be removed by finishing procedures. For intra-oral application, the material is for transitory use, that is, it can be applied in the mouth only during the molding and sculpture process.

Warnings, precautions and restrictions: Professional use. Handle the product wearing gloves. The product must not be used in patients who show allergic reactions to any of its components. To prevent the product from polymerizing, close the syringe immediately after use. The packaging should be protected from heat and direct exposure to sunlight. The Material Safety Data Sheet for this product can be obtained by telephone. Disposal: Waste must be disposed of in an industrial or sanitary landfill, in compliance with local legislation in force. Do not dispose of waste in sewers, waterways or household waste bin. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. * May cause allergic skin reactions. * Toxic to aquatic organisms with long lasting effects. Contaminated clothing must not leave from its workplace. Avoid release to the environment. In case of contact with skin: Wash with plenty of soap and water. Specific treatment. In case of skin irritation or rash: Get medical advice / attention. Remove all contaminated clothing and wash before reuse. Put away the spilled material. Dispose of contents / container in accordance with local regulations. **Storage:** To be stored in a cool dry place, protected from sunlight. Do not store above 30°C.

Composition: Methacrylate and activators.

Ref.	Código/Code	Conteúdo/Contenido/ Content:
0.1261	0012600006	3g pasta/paste

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