EN SeptoCure Instrumentos to condense composites

SeptoCure is indicated to condense composite resins in clinical procedures, since it allows to light cure through their tips while the composites are being hold under pressure on the bottom of the cavity during polymerization, reducing the risk of marginal leakage. Two instruments with four different tips, developed specifically for the application of this technique.

Use: Select among the four available tips the one which best adapts to the proximal box. Take the first increment of resin to the cavity, keeping it under pressure while being polymerized. Repeat the procedure, layer by layer, as often as necessary.

Warnings, precautions and restrictions: Professional use. If there is proven allergy to any of its components, the product must not be used. The instrument should be sterilized in autoclave before use. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amilase, Protease, Lipase and Carbohydrase) and with pure pH puro between 6.5 - 7.5. Prepare the dilution solution according to the manufacturer's instructions. The water should be potable. For best effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the products as follows: immerse them in contact with the solution for 5 minutes. Remove parts and rinse under running water for at least 10 seconds. Dry the piece with absorbent paper or disposable wipes. **Sterilization:** The instrument should be autoclaved (126°C, 16 min.), after each patient. Sterilization performed with any method other than autoclave can damage the product. **Reprocessing Limits:** This is a reusable device. The manufacturer guarantees up to 300 reuses.

Inspection: Check the product before each use, observing whether there is damage to the tips or cracks along the instrument, and do not use in such cases. Cracks are indicative of material fragility, which can cause the product to break.

Composition: Polysulfone.
Storage: To be stored in a cool and dry place, protected from sunlight.
Disposal: Dispose in accordance with applicable local requirements.
Presentación: 0105010002 - SeptoCure – 2 instruments.

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